

## Objective

To test performance and calibrate the instrument for vFC™ counting and sizing. For Protocol 0 on an instrument that is appropriately configured, we would expect to see approximately 50,000 events with a diameter distribution that ranges between ~50 – 300 nm, and a median of 100-130 nm and mean of ~120-150 nm.

## Materials

- a. Gloves
- b. Staining Buffer, 8 mL
- c. Lipo100™ Standard (10x)
- d. vFRed™ membrane stain (100x)
- e. Microwell plate

## Procedure

1. Wear gloves for all manipulation of samples and reagents.

### Prepare Working Solution

2. Prepare 10x vFRed™ working solution by adding 8 µL vFRed (100x) to 72 µL Staining Buffer.

### Prepare Samples

3. Prepare EV samples (See Protocol 0 Plate Map below):
  - a. Pipet 40 µL of Staining Buffer into wells A1-H1
  - b. Add 5 µL Lipo100™ (10x) to wells C1, D1, G1, and H1.
  - c. Add 5µL of staining buffer into wells A1, B1, E1, and F1
4. Add 5 µL 10x vFRed™ to all wells except A1 and E1.
5. Mix well, incubate 1 hour at RT in the dark.

| Table 1. Staining reactions |              |              |
|-----------------------------|--------------|--------------|
| Reagent                     | CytoFlex     | CellStream   |
| Buffer                      | 40 uL        | 35 uL        |
| Lipo100                     | 5 uL         | 10 uL        |
| vFRed                       | 5 uL         | 5 uL         |
| <b>Total</b>                | <b>50 uL</b> | <b>50 uL</b> |

### Dilute and Read

6. Dilute the staining reaction according to Table 2.

For the CytoFlex®

- a. Add 145ul of staining buffer to row 2 and 291uL of staining buffer to row 3.
- b. Transfer 5 µL of stained sample (row 1) into a well containing 145 µL of Staining Buffer, mix well by pipetting up and down (this will be dilution 1). Be careful to avoid foaming.
- c. Transfer 9 µL of Dilution 1 into a well containing 291 µL of Staining Buffer, mix well (this will be dilution 2).
- d. Run Dilution 2 on CytoFlex for fixed time (120 seconds) at fixed flow rate (High, 60 µL/min)

For the CellStream®

- e. Add 298.5ul of staining buffer to row 2.
- f. Transfer 1.5 µL of stained sample (row 1) into a well containing 298.5 µL of Staining Buffer, mix well by pipetting up and down. Be careful to avoid foaming.
- g. Run on CellStream for fixed time (120 seconds) at fixed flow rate (Low, 3.66 µL/min)

| Table 2. Post-stain dilution and run |                   |                    |
|--------------------------------------|-------------------|--------------------|
| Reagent                              | CytoFlex          | CellStream         |
| Staining                             | 50 uL             | 50 uL              |
| Dilution 1                           | 5 µL →<br>145 uL  | 1.5 µL →<br>300 uL |
| <u>Dilution 2</u>                    | 9 µL →<br>291 uL  |                    |
| <b>Dilution factor</b>               | 1000              | 200                |
| <b>Run speed</b>                     | High<br>60 uL/min | Slow<br>3.7 uL/min |

## Data Analysis

- Open the vFC™ Analysis layout with FCS Express Reader.  
**Note: The vFC Analysis Layout and Chapter 4 below have additional notes and tips to guide the data analysis.**
- From the Data List, click the Add File (+) and navigate to the data directory and select the Protocol 0 data files).
  - Select a **Buffer + vFRed™** data file and click Change Data On All Plots. Inspect the Gating Plots. Adjust Time Gate, Pulse Gate, and Vesicle Gate as needed to minimize low- and high-scatter backgrounds.
  - Select a **Lipo100™** data file and click Change Data On All Plots. Inspect the Gating Plots. Adjust Gates as needed to select vFRed™-positive events.
- Run the Vesicle Calibration Tool or Copy the Lipo100 vFRed™ intensity data to the clipboard, and paste into the vFC Size Calibration Tool located under the *Data Analysis* tab at [www.cellarcus.com/products/vfc-base-kit](http://www.cellarcus.com/products/vfc-base-kit). Enter the parameters from calibration results into the Surface Area parameter via Tools>Transforms>Parameter Math.
- Inspect the resulting size distribution. The Lipo100 vesicle estimated size should range between ~50 – 300 nm, and a median 100-130 nm and mean of ~120-150 nm.
- Export plots and statistics via Batch Processing.

## Protocol 0 Plate Map

|   | Staining wells         | Dilute 1<br>Don't run | Dilute 2<br>Run on FC | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|------------------------|-----------------------|-----------------------|---|---|---|---|---|---|----|----|----|
| A | 1<br>Buffer only       | 2<br>5 uL into 145 uL | 3<br>9 uL into 291 uL |   |   |   |   |   |   |    |    |    |
| B | Buffer + vFRed         | 5 uL into 145 uL      | 9 uL into 291 uL      |   |   |   |   |   |   |    |    |    |
| C | Lipo100™ Size Standard | 5 uL into 145 uL      | 9 uL into 291 uL      |   |   |   |   |   |   |    |    |    |
| D | Lipo100™ Size Standard | 5 uL into 145 uL      | 9 uL into 291 uL      |   |   |   |   |   |   |    |    |    |
| E | Buffer only            | 5 uL into 145 uL      | 9 uL into 291 uL      |   |   |   |   |   |   |    |    |    |
| F | Buffer + vFRed         | 5 uL into 145 uL      | 9 uL into 291 uL      |   |   |   |   |   |   |    |    |    |
| G | Lipo100™ Size Standard | 5 uL into 145 uL      | 9 uL into 291 uL      |   |   |   |   |   |   |    |    |    |
| H | Lipo100™ Size Standard | 5 uL into 145 uL      | 9 uL into 291 uL      |   |   |   |   |   |   |    |    |    |

Developed, Manufactured, and Distributed By

**Cellarcus Biosciences, Inc.**

Telephone: +1 (858) 239-2100

Customer Care: [cellarcus@cellarcus.com](mailto:cellarcus@cellarcus.com)

Technical Support: [technical@cellarcus.com](mailto:technical@cellarcus.com)

Notes